

Nur für Forschungszwecke

TCEAL3 Polyklonaler Antikörper

Katalog-Nr.:10732-2-AP



Allgemeine Informationen

Katalog-Nr.:
10732-2-AP

Größe:

150ul, Konzentration: 147 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG1039

GenBank-Zugangsnummer:

BC007622

GeneID (NCBI):

85012

Vollständiger Name:

transcription elongation factor A (SII)-
like 3

Berechnete Masse:

23 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

IHC,ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human ovary tumor using 10732-2-AP (TCEAL3 antibody) at dilution of 1:100 (under 40x lens).