

Nur für Forschungszwecke

CD58 Polyklonaler Antikörper

Katalog-Nr.:10878-1-AP



Allgemeine Informationen

Katalog-Nr.:
10878-1-AP

Größe:
150ul , Konzentration: 500 µg/ml von
Nanodrop und 307 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG1165

GenBank-Zugangsnummer:
BC005930

GeneID (NCBI):
965

Vollständiger Name:
CD58 molecule

Berechnete Masse:
28 kDa

Beobachtete Masse:
55-70 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB : Raji-Zellen, HeLa-Zellen, humanes
Lymphgewebe, Jurkat-Zellen, TF-1-Zellen

Hintergrundinformationen

CD58 a heavily glycosylated protein of 55-70 kDa, is expressed on a broad range of hematopoietic and non-hematopoietic cells. CD58 is involved in immune recognition of tumor cells via binding of the CD2 receptor expressed on cytotoxic T cells. The protein is localized to the plasma membrane.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various cells were subjected to SDS PAGE followed by western blot with 10878-1-AP (CD58 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.