

Nur für Forschungszwecke

ABTB1 Polyklonaler Antikörper

Katalog-Nr.:11095-1-AP



Allgemeine Informationen

Katalog-Nr.: 11095-1-AP	GenBank-Zugangsnummer: BC011858	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 240 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 80325	
Wirt: Kaninchen	Vollständiger Name: ankyrin repeat and BTB (POZ) domain containing 1	
Isotyp: IgG	Berechnete Masse: 54 kDa	
Immunogen Katalognummer: AG1574		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

Proteins with a bric-à-brac, tramtrack, broad-complex/Poxvirus zinc fingers (BTB/POZ) domain are involved in a broad variety of biological processes, including DNA binding, organization of macromolecular structures and regulation of gene transcription [PMID:11959095]. BPOZ contains an N-terminal ankyrin repeat, a bipartite nuclear localization signal and two BTB/POZ domains. It may involve in the development processes [PMID:10891360]. ABTB1 exists various isoforms and range of molecular weight of isoforms is 45-54 kDa.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

