

Nur für Forschungszwecke

PITX2 Polyklonaler Antikörper

Katalog-Nr.:11106-1-AP



Allgemeine Informationen

Katalog-Nr.:
11106-1-AP

Größe:
150ul , Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG1583

GenBank-Zugangsnummer:
BC013998

GeneID (NCBI):
5308

Vollständiger Name:
paired-like homeodomain 2

Berechnete Masse:
36 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

PITX2 belongs to the bicoid class of homeodomain transcription factors, which play essential roles in embryonic development. Various PITX2 isoforms are involved in a wide variety of developmental programs, including left-right signaling. It controls cell proliferation in a tissue-specific manner and is involved in morphogenesis. During embryonic development, PITX2 involves in the expansion of muscle progenitors. In addition, it may play a role in the proper localization of asymmetric organs such as the heart and stomach.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

