

Nur für Forschungszwecke

EIF2AK4 Polyklonaler Antikörper

Katalog-Nr.: 11174-1-AP



Allgemeine Informationen

Katalog-Nr.:

11174-1-AP

Größe:

150ul, Konzentration: 127 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG1654

GenBank-Zugangsnummer:

BC009350

GeneID (NCBI):

440275

Vollständiger Name:

eukaryotic translation initiation
factor 2 alpha kinase 4

Berechnete Masse:

187 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

IHC, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human prostate cancer using 11174-1-AP (EIF2AK4 antibody) at dilution of 1:50 (under 10x lens).