

PRMT2 Polyklonaler Antikörper

Katalog-Nr.: 11269-1-AP

Allgemeine Informationen

Katalog-Nr.:
11269-1-AP

Größe:
150ul, Konzentration: 140 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG1806

GenBank-Zugangsnummer:
BC000727

GeneID (NCBI):
3275

Vollständiger Name:
protein arginine methyltransferase 2

Berechnete Masse:
49 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Positivkontrollen:

WB: humanes Plazenta-Gewebe,

Hintergrundinformationen

PRMT2 (also known as HRMT1L1) is a member of the arginine methyltransferase family which catalyzes arginine methylation and regulates diverse cellular processes, including transcription, translation, DNA repair, and RNA processing. Recently PRMT2 has been identified as a pro-tumorigenic factor in malignant gliomas.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

