

GDE1 Polyklonaler Antikörper

Katalog-Nr.: 11272-2-AP

1 Publikationen

Allgemeine Informationen

Katalog-Nr.:
11272-2-AP

Größe:
150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG1821

GenBank-Zugangsnummer:
BC014981

GeneID (NCBI):
51573

Vollständiger Name:
glycerophosphodiester
phosphodiesterase 1

Berechnete Masse:
38 kDa

Beobachtete Masse:
55-66 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

In Publikationen genannte Anwendungen:
WB

Getestete Reaktivität:
Human, Maus, Ratte

Zitierte Arten:
Human

Hintergrundinformationen

GDE1(glycerophosphodiester phosphodiesterase 1) is also named as MIR16 and belongs to the glycerophosphoryl diester phosphodiesterase family. It is an 43 kDa integral membrane glycoprotein, which remains associated with membranes after alkaline (pH 11) extraction and is sensitive to endoglycosidase H and most likely involved in the remodeling of membrane PIs and regulation of membrane phospholipid composition, and hence in intracellular membrane trafficking, as well as in the modulation of G proteinmediated signaling(PMID:12576545). The full length protein has two glycosylation sites.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Sasha M Milbeck	32602905	Invest Ophthalmol Vis Sci	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

