

Nur für Forschungszwecke

# Livin Polyklonaler Antikörper

Katalog-Nr.:11294-2-AP



## Allgemeine Informationen

Katalog-Nr.:  
11294-2-AP

Größe:  
150ul , Konzentration: 133 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG1861

GenBank-Zugangsnummer:  
BC014475

GeneID (NCBI):  
79444

Vollständiger Name:  
baculoviral IAP repeat-containing 7

Berechnete Masse:  
33 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:  
IHC,ELISA

Getestete Reaktivität:  
Human

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil  
Lagerungspuffer:  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human colon cancer using 11294-2-AP (Livin antibody) at dilution of 1:50 (under 10x lens).