

Nur für Forschungszwecke

SNAPC1 Polyklonaler Antikörper

Katalog-Nr.: 11297-1-AP



Allgemeine Informationen

Katalog-Nr.:
11297-1-AP

Größe:
150ul, Konzentration: 233 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG1812

GenBank-Zugangsnummer:
BC019038

GeneID (NCBI):
6617

Vollständiger Name:
small nuclear RNA activating
complex, polypeptide 1, 43kDa

Berechnete Masse:
43 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
IHC, ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human colon cancer using 11297-1-AP (SNAPC1 antibody) at dilution of 1:50 (under 10x lens).