

Nur für Forschungszwecke

# DEF6 Polyklonaler Antikörper

Katalog-Nr.:11369-1-AP

Vorgestelltes Produkt



## Allgemeine Informationen

Katalog-Nr.:  
11369-1-AP

Größe:

150ul, Konzentration: 350 µg/ml von  
Nanodrop und 200 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:  
AG1925

GenBank-Zugangsnummer:

BC017504

GeneID (NCBI):

50619

Vollständiger Name:

differentially expressed in FDCP 6  
homolog (mouse)

Berechnete Masse:

74 kDa

Beobachtete Masse:

74 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

IP 0.5-4.0 µg für IP und 1:500-1:1000  
für WB

IHC 1:50-1:500

## Anwendungen

Geprüfte Anwendungen:

IHC, IP, WB, ELISA

Getestete Reaktivität:

Human, Maus

**Hinweis-IHC: Antigendemaskierung mit TE-  
Puffer pH 9,0 empfohlen. (\*) Wahlweise  
kann die Antigendemaskierung auch mit  
Citratpuffer pH 6,0 erfolgen.**

Positivkontrollen:

WB: Jurkat-Zellen,

IP: Jurkat-Zellen,

IHC: Mausmilzgewebe, humanes Tonsillitisgewebe

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.


## Ausgewählte Validierungsdaten




Jurkat cells were subjected to SDS PAGE followed by western blot with 11369-1-AP (DEF6 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.




WB result of DEF6 antibody (11369-1-AP; 1:4000; incubated at room temperature for 1.5 hours) with sh-Control and sh-DEF6 transfected Jurkat cells.



Immunohistochemical analysis of paraffin-embedded mouse spleen tissue slide using 11369-1-AP (DEF6 antibody) at dilution of 1:200 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse spleen tissue slide using 11369-1-AP (DEF6 antibody) at dilution of 1:200 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



IP Result of anti-DEF6 (IP:11369-1-AP, 4ug; Detection:11369-1-AP 1:600) with Jurkat cells lysate 2400ug.