

Nur für Forschungszwecke

FTSJ1 Polyklonaler Antikörper

Katalog-Nr.:11620-1-AP



Allgemeine Informationen

Katalog-Nr.: 11620-1-AP	GenBank-Zugangsnummer: BC023584	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 900 µg/ml von Nanodrop;	GeneID (NCBI): 24140	Empfohlene Verdünnungen: WB 1:1000-1:4000 IHC 1:500-1:2000 IF 1:200-1:800
Wirt: Kaninchen	Vollständiger Name: FtsJ homolog 1 (E. coli)	
Isotyp: IgG	Berechnete Masse: 329 aa, 36 kDa	
Immunogen Katalognummer: AG2199	Beobachtete Masse: 36 kDa	

Anwendungen

Geprüfte Anwendungen:
IF, IHC, WB, ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hinweis-IHC: Antigenmaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigenmaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Positivkontrollen:

WB: HeLa-Zellen, HepG2-Zellen, U-87 MG-Zellen

IHC: humanes Lebergewebe, Maushirngewebe

IF: A431-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




HeLa cells were subjected to SDS PAGE followed by western blot with 11620-1-AP (FTSJ1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed A431 cells using FTSJ1 antibody (11620-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 11620-1-AP (FTSJ1 antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 11620-1-AP (FTSJ1 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).