

Nur für Forschungszwecke

MCF2L Polyklonaler Antikörper

Katalog-Nr.: 11659-1-AP



Allgemeine Informationen

Katalog-Nr.:
11659-1-AP

Größe:

150ul, Konzentration: 93 µg/ml durch
die Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG2223

GenBank-Zugangsnummer:

BC020208

GeneID (NCBI):

23263

Vollständiger Name:

MCF.2 cell line derived transforming
sequence-like

Berechnete Masse:

984aa, 112 kDa; 1108aa, 124 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

IHC, ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human gliomas using 11659-1-AP (MCF2L antibody) at dilution of 1:100 (under 10x lens).