

Nur für Forschungszwecke

NFATC2IP Polyklonaler Antikörper

Katalog-Nr.:11667-1-AP



Allgemeine Informationen

Katalog-Nr.: 11667-1-AP	GenBank-Zugangsnummer: BC021551	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 247 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 84901	
Wirt: Kaninchen	Vollständiger Name: nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 2 interacting protein	
Isotyp: IgG	Berechnete Masse: 138 aa, 15 kDa	
Immunogen Katalognummer: AG2217	Beobachtete Masse: 50 kDa	

Anwendungen

Geprüfte Anwendungen:
IHC, WB,ELISA

Getestete Reaktivität:
Human, Maus

Hintergrundinformationen

NFATC2IP, also named as NIP45, is a 419 amino acid protein, which contains one ubiquitin-like domain. TRAF1 is associated with a fraction of NFATC2IP in the cytoplasm and prevents its translocation to the nucleus. NFATC2IP in T-helper 2 (Th2) cells, regulates the magnitude of NFAT-driven transcription of a specific subset of cytokine genes, including IL3, IL4, IL5 and IL13, but not IL2.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 11667-1-AP (NFATC2IP Antibody) at dilution of 1:100 incubated at 4 degree celsius over night.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 11667-1-AP (NFATC2IP Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 11667-1-AP (NFATC2IP Antibody) at dilution of 1:50 (under 40x lens).