

Nur für Forschungszwecke

# SFRS6 Polyklonaler Antikörper

Katalog-Nr.: 11772-1-AP

3 Publikationen



## Allgemeine Informationen

<b>Katalog-Nr.:</b> 11772-1-AP	<b>GenBank-Zugangsnummer:</b> BC006832	<b>Reinigungsmethode:</b> Antigen-Affinitätsreinigung
<b>Größe:</b> 150ul, Konzentration: 450 µg/ml von Nanodrop;	<b>GeneID (NCBI):</b> 6431	<b>Empfohlene Verdünnungen:</b> WB 1:500-1:3000
<b>Wirt:</b> Kaninchen	<b>Vollständiger Name:</b> splicing factor, arginine/serine-rich 6	
<b>Isotyp:</b> IgG	<b>Berechnete Masse:</b> 344 aa, 40 kDa	
<b>Immunogen Katalognummer:</b> AG2412	<b>Beobachtete Masse:</b> 40 kDa	

## Anwendungen

<b>Geprüfte Anwendungen:</b> WB, ELISA	<b>Positivkontrollen:</b> WB : HeLa-Zellen, MCF-7-Zellen
<b>In Publikationen genannte Anwendungen:</b> WB	
<b>Getestete Reaktivität:</b> Human	
<b>Zitierte Arten:</b> Human, Maus	

## Hintergrundinformationen

SR proteins have a characteristic C-terminal Ser/Arg-rich repeat (RS domain) of variable length and constitute a family of highly conserved nuclear phosphoproteins that can function as both essential and alternative pre-mRNA splicing factors [PMID:7556075]. SRF56 has a pivotal role in the TF exon 5 definition, an event essential for flTF biosynthesis. Besides, up- or down regulation of SRSF6 alters the splicing of several tumor suppressors and oncogenes to generate the oncogenic isoforms and reduce the tumor suppressive isoforms [PMID:23132731].

## Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Norihiro Okada	34082065	Gene	WB
Yuhan Chen	29942010	Nat Commun	WB
Norihiro Okada	35515992	Data Brief	WB

## Lagerung


**Lagerungsbedingungen:**  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil  
**Lagerungspuffer:**  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 11772-1-AP (SFRS6 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.