

Nur für Forschungszwecke

RSAD2 Polyklonaler Antikörper

Katalog-Nr.: 11833-1-AP

12 Publikationen



Allgemeine Informationen

Katalog-Nr.:
11833-1-AP

Größe:
150ul, Konzentration: 440 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG2409

GenBank-Zugangsnummer:
BC017969

GeneID (NCBI):
91543

Vollständiger Name:
radical S-adenosyl methionine
domain containing 2

Berechnete Masse:
361 aa, 42 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

In Publikationen genannte Anwendungen:
WB

Getestete Reaktivität:
Human, Maus, Ratte

Zitierte Arten:
Human, Maus, Rind

Hintergrundinformationen

RSAD2 (radical S-adenosyl methionine domain-containing protein 2), also known as CIG5 (cytomegalovirus-induced gene 5 protein), vlg1, viperin or CIG33, displays antiviral effect against HIV-1 virus, hepatitis C virus, human cytomegalovirus, and aphaviruses. Expression of the protein can be induced by interferon. Western blot analysis of IFNG-treated cells revealed 43 kD/35kd bands.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Timothy J Grunkemeyer	34029588	J Biol Chem	WB
Arti B Dumbrepatil	30872404	J Biol Chem	
Ayesha M Patel	35831548	Sci Rep	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

