

Nur für Forschungszwecke

Langerin Polyklonaler Antikörper

Katalog-Nr.:11841-1-AP



Allgemeine Informationen

Katalog-Nr.: 11841-1-AP	GenBank-Zugangsnummer: BC022278	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 180 µg/ml von Nanodrop und 133 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 50489	Empfohlene Verdünnungen: WB 1:500-1:1000 IHC 1:3000-1:8000
Wirt: Kaninchen	Vollständiger Name: CD207 molecule, langerin	
Isotyp: IgG	Berechnete Masse: 328 aa, 37 kDa	
Immunogen Katalognummer: AG2479	Beobachtete Masse: 40 kDa	

Anwendungen

Geprüfte Anwendungen: IHC, WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Human, Maus	WB: humanes Hautgewebe, Maushautgewebe
Hinweis-IHC: Antigenmaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigenmaskierung auch mit Citratpuffer pH 6,0 erfolgen.	IHC: humanes Hautkrebsgewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten



human skin tissue were subjected to SDS PAGE followed by western blot with 11841-1-AP (Langerin antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human skin cancer tissue slide using 11841-1-AP (Langerin antibody) at dilution of 1:6000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human skin cancer tissue slide using 11841-1-AP (Langerin antibody) at dilution of 1:6000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).