

Nur für Forschungszwecke

VPREB3 Polyklonaler Antikörper

Katalog-Nr.: 11915-1-AP



Allgemeine Informationen

Katalog-Nr.: 11915-1-AP	GenBank-Zugangsnummer: BC020666	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 133 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 29802	
Wirt: Kaninchen	Vollständiger Name: pre-B lymphocyte 3	
Isotyp: IgG	Berechnete Masse: 123 aa, 14 kDa	
Immunogen Katalognummer: AG2528		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

VPREB3 (formerly named 8HS20) is a protein that associates with mu chains in the endoplasmic reticulum of pre-B cell lines (PMID: 8491176; 22684248). VPREB3 is normally expressed by precursor B cells in bone marrow and by a subset of normal germinal center B cells in secondary lymphoid organs. It has been reported that VPREB3 expression is associated with B-cell tumors with c-MYC abnormalities (PMID: 20823132).

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

