

Allgemeine Informationen

Katalog-Nr.:

11930-1-AP

Größe:

150ul, Konzentration: 900 µg/ml von Nanodrop und 453 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG2464

GenBank-Zugangsnummer:

BC022854

GeneID (NCBI):

64231

Vollständiger Name:

membrane-spanning 4-domains, subfamily A, member 6A

Berechnete Masse:

248 aa, 27 kDa

Beobachtete Masse:

40-45 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:3000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : K-562-Zellen, humanes Hirngewebe, Jurkat-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



K-562 cells were subjected to SDS PAGE followed by western blot with 11930-1-AP (MS4A6A antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.