

Nur für Forschungszwecke

PTMS Polyklonaler Antikörper

Katalog-Nr.: 12106-1-AP



Allgemeine Informationen

Katalog-Nr.: 12106-1-AP	GenBank-Zugangsnummer: BC017025	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 133 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 5763	Empfohlene Verdünnungen: WB 1:200-1:1000 IF 1:10-1:100
Wirt: Kaninchen	Vollständiger Name: parathymosin	
Isotyp: IgG	Berechnete Masse: 12 kDa	
Immunogen Katalognummer: AG2751	Beobachtete Masse: 10-12 kDa	

Anwendungen

Geprüfte Anwendungen: IF, WB, ELISA	Positivkontrollen: WB : U-937-Zellen, IF : A549-Zellen,
Getestete Reaktivität: Human	

Hintergrundinformationen

PTMS, also named as Parathymosin, is a small nuclear protein that can physically interact with glucocorticoid receptors (GR) and histone H1, respectively, and can serve as a coactivator of GR (PMID: 16150697; PMID: 15716277). It is widely distributed in mammalian tissues. Catalog# 12106-1-AP is a rabbit polyclonal antibody raised against the full-length of human PTMS. The MW of this protein is 12 kDa, and this antibody specially recognises the 12 kDa protein.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



U-937 cells were subjected to SDS PAGE followed by western blot with 12106-1-AP (PTMS antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of A549 cells using 12106-1-AP (PTMS antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).