

Nur für Forschungszwecke

VPS33B Polyklonaler Antikörper

Katalog-Nr.:12195-1-AP

Vorgestelltes Produkt

11 Publikationen



Allgemeine Informationen

Katalog-Nr.: 12195-1-AP	GenBank-Zugangsnummer: BC016445	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 300 µg/ml von Nanodrop;	GeneID (NCBI): 26276	Empfohlene Verdünnungen: WB 1:500-1:2000 IHC 1:50-1:500 IF 1:10-1:100
Wirt: Kaninchen	Vollständiger Name: vacuolar protein sorting 33 homolog B (yeast)	
Isotyp: IgG	Berechnete Masse: 617 aa, 71 kDa	
Immunogen Katalognummer: AG2833	Beobachtete Masse: 71 kDa	

Anwendungen

Geprüfte Anwendungen:

IF, IHC, WB, ELISA

In Publikationen genannte Anwendungen:

ChIP, CoIP, IF, IHC, IP, WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human, Maus, Zebrafisch

Hinweis-IHC: Antigenmaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigenmaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Positivkontrollen:

WB: HEK-293-Zellen, HeLa-Zellen, L02-Zellen, Maushodengewebe, Maus-Skelettmuskelgewebe

IHC: humanes Kolonkarzinomgewebe, humanes Eierstockgewebe, humanes Hirngewebe, humanes Hodengewebe, humanes Lungengewebe, humanes Medulloblastom-Gewebe, humanes Milzgewebe, humanes Nierengewebe, humanes Plazenta-Gewebe

IF: MCF-7-Zellen,

Hintergrundinformationen

VPS33B, a homolog of yeast class C vacuolar protein sorting (vps) protein Vps33p, belongs to the STXBP/unc-18/SEC1 family. It may play a role in vesicle-mediated protein trafficking to lysosomal compartments and in membrane docking/fusion reactions of late endosomes/lysosomes. VPS33B mediates phagolysosomal fusion in macrophages (PMID: 18474358). Defects in VPS33B account for most cases of arthrogryposis-renal dysfunction-cholestasis (ARC) syndrome, which is a multisystem disorder associated with abnormalities in polarized liver and kidney cells (PMID: 20190753).

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Rik van der Kant	26463206	J Biol Chem	WB, IF
Yiyu Chen	31125123	Int J Cancer	WB, IHC
Yu Liu	32385361	Sci Rep	IF

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




HEK-293 cells were subjected to SDS PAGE followed by western blot with 12195-1-AP (VPS33B Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of MCF-7 cells, using VPS33B antibody 12195-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 12195-1-AP (VPS33B Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 12195-1-AP (VPS33B Antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).