

Nur für Forschungszwecke

INPP4A Polyklonaler Antikörper

Katalog-Nr.: 12665-1-AP



Allgemeine Informationen

Katalog-Nr.:
12665-1-AP

Größe:
150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG3346

GenBank-Zugangsnummer:
BC028361

GeneID (NCBI):
3631

Vollständiger Name:
inositol polyphosphate-4-
phosphatase, type I, 107kDa

Berechnete Masse:
972 aa, 109 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

IP 0.5-4.0 µg für IP und 1:200-1:1000
für WB

Anwendungen

Geprüfte Anwendungen:
IP, ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Positivkontrollen:

IP: Maushirngewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



IP result of anti-INPP4A (IP:12665-1-AP, 4ug;
Detection:12665-1-AP 1:300) with mouse brain
tissue lysate 3200 ug.