

Nur für Forschungszwecke

VASH1 Polyklonaler Antikörper

Katalog-Nr.:12730-1-AP

Vorgestelltes Produkt

2 Publikationen



Allgemeine Informationen

Katalog-Nr.:
12730-1-AP

Größe:
150ul, Konzentration: 350 µg/ml von
Nanodrop und 287 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG3435

GenBank-Zugangsnummer:
BC009031

GeneID (NCBI):
22846

Vollständiger Name:
vasohibin 1

Berechnete Masse:
365 aa, 41 kDa

Beobachtete Masse:
36 and 41 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human

Zitierte Arten:

Zebrafisch

Positivkontrollen:

WB : humanes Hirngewebe, HeLa-Zellen, MCF-7-Zellen

Hintergrundinformationen

VASH1, also named as vasohibin 1, is a angiogenesis inhibitor induced by the stimulation with VEGF. High expression of VASH1 has been researched in some tumors as a result of VEGF induction. It has been suggested that VASH-1 is clinically significant in bladder cancer, colorectal cancer and NSCLC. VASH-1 may be a potential marker of gastric cancer and a target for vascular targeting therapy. 12730-1-AP detects 41 kDa form secreted protein and 36 kDa form protein accumulated within the cells or pericellular milieu in SDS-PAGE.(PMID:19057892, 29042936, 29552154, 26666821)

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Marta Bastos de Oliveira	33547133	Development	WB
Ana C Almeida	35385739	Cell Rep	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



human brain tissue were subjected to SDS PAGE followed by western blot with 12730-1-AP (VASH1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.