

Nur für Forschungszwecke

# RIC3 Polyklonaler Antikörper

Katalog-Nr.: 12779-1-AP



## Allgemeine Informationen

Katalog-Nr.:

12779-1-AP

Größe:

150ul, Konzentration: 120 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG3479

GenBank-Zugangsnummer:

BC022455

GeneID (NCBI):

79608

Vollständiger Name:

resistance to inhibitors of  
cholinesterase 3 homolog (C.  
elegans)

Berechnete Masse:

368 aa, 41 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:

IHC, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human prostate cancer using 12779-1-AP (RIC3 antibody) at dilution of 1:100 (under 10x lens).