

Nur für Forschungszwecke

TSPAN3 Polyklonaler Antikörper

Katalog-Nr.:13036-1-AP



Allgemeine Informationen

Katalog-Nr.: 13036-1-AP	GenBank-Zugangsnummer: BC011206	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 133 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 10099	
Wirt: Kaninchen	Vollständiger Name: tetraspanin 3	
Isotyp: IgG	Berechnete Masse: 253 aa, 28 kDa	
Immunogen Katalognummer: AG3689		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

Tetraspanin-3 (TSPAN3, also known as TM4SF8) is a member of the tetraspanin (TM4SF) family. Tetraspanins, so named because they contain four transmembrane domains, are frequently expressed at the cell surface in association with each other and with other molecules, and they function to regulate cell adhesion, migration, proliferation, and differentiation. TSPAN3 regulates the proliferation and migration of oligodendrocytes, a process essential for normal myelination and repair.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

