

Nur für Forschungszwecke

MPZL1 (BC019890, 147aa) Polyklonaler Antikörper



Katalog-Nr.: 13158-1-AP

1 Publikationen

Allgemeine Informationen

Katalog-Nr.:

13158-1-AP

Größe:

150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG3791

GenBank-Zugangsnummer:

BC019890

GeneID (NCBI):

9019

Vollständiger Name:

myelin protein zero-like 1

Berechnete Masse:

147 aa, 16 kDa

Beobachtete Masse:

16 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human, Maus, Ratte

Zitierte Arten:

Human

Hintergrundinformationen

Myelin protein zero-like 1 (MPZL1), also known as Protein zero-related, is a transmembrane glycoprotein involved in extracellular matrix-induced signal transduction (PMID: 9792637; 11751924). Its extracellular segment has significant sequence homology to myelin P0, whereas its intracellular portion has two tyrosine phosphorylation sites resembling ITIMs (PMID: 9792637). Alternative splicing of the MPZL1 gene results in multiple transcript variants encoding different isoforms. This antibody raised against MPZL1 (BC019890, 147aa).

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Yunyao Ye	32308413	Onco Targets Ther	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

