

Nur für Forschungszwecke

ETV6 Polyklonaler Antikörper

Katalog-Nr.:13182-1-AP



Allgemeine Informationen

Katalog-Nr.:

13182-1-AP

Größe:

150ul , Konzentration: 133 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG3852

GenBank-Zugangsnummer:

BC043399

GeneID (NCBI):

2120

Vollständiger Name:

ets variant 6

Berechnete Masse:

452 aa, 53 kDa

Beobachtete Masse:

57 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:4000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : A431-Zellen, HeLa-Zellen, Jurkat-Zellen, Raji-Zellen

Hintergrundinformationen

ETV6, the E26 transformation-specific (ETS) family transcription repressor and tumor suppressor variant 6 gene, resides on chromosome 12p13, encodes a 57 kDa ETS family transcription factor with 3 highly conserved functional domains(PMID: 32841218). ETV6 is a well-studied hematopoietic transcription factor reported to play a key functional role in hematopoiesis, megakaryopoiesis, cell growth, and differentiation.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 13182-1-AP (ETV6 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.