

Nur für Forschungszwecke

PAX9 Polyklonaler Antikörper

Katalog-Nr.: 13195-1-AP



Allgemeine Informationen

Katalog-Nr.:

13195-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG3897

GenBank-Zugangsnummer:

BC001159

GeneID (NCBI):

5083

Vollständiger Name:

paired box 9

Berechnete Masse:

341 aa, 36 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

PAX9, a member of the paired box-containing gene family, is closely related in its paired domain to PAX1. PAX9 is essential for the development of a variety of organs and skeletal elements, particularly the initial stage of odontogenesis. Mutations in PAX9 can cause selective tooth agenesis in humans. PAX9 is reported as a sensitive marker for deregulated differentiation of oesophageal keratinocytes, and it has a role in the normal differentiation process of internal stratified squamous epithelia.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

