

Nur für Forschungszwecke

ATCAY Polyklonaler Antikörper

Katalog-Nr.: 13257-2-AP



Allgemeine Informationen

Katalog-Nr.:
13257-2-AP

Größe:
150ul, Konzentration: 280 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG3658

GenBank-Zugangsnummer:
BC026217

GeneID (NCBI):
85300

Vollständiger Name:
ataxia, cerebellar, Cayman type

Berechnete Masse:
371 aa, 42 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse brain tissue were subjected to SDS PAGE followed by western blot with 13257-2-AP (ATCAY antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.