

Nur für Forschungszwecke

# LENG8 Polyklonaler Antikörper

Katalog-Nr.: 13336-1-AP

1 Publikationen



## Allgemeine Informationen

Katalog-Nr.:  
13336-1-AP

Größe:  
150ul, Konzentration: 507 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG4153

GenBank-Zugangsnummer:  
BC028048

GeneID (NCBI):  
114823

Vollständiger Name:  
leukocyte receptor cluster (LRC)  
member 8

Berechnete Masse:  
800 aa, 88 kDa

Beobachtete Masse:  
86-88 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:  
WB 1:500-1:3000

## Anwendungen

Geprüfte Anwendungen:  
FC, WB, ELISA

In Publikationen genannte Anwendungen:  
WB

Getestete Reaktivität:  
Human, Maus, Ratte

Zitierte Arten:  
Human

Positivkontrollen:

WB : A2780-Zellen, Jurkat-Zellen

## Hintergrundinformationen

### Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Xiaojuan Ye	26884854	Int J Clin Exp Pathol	WB

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Ausgewählte Validierungsdaten



A2780 cells were subjected to SDS PAGE followed by western blot with 13336-1-AP (LENG8 antibody) at dilution of 1:1200 incubated at room temperature for 1.5 hours.



1X10<sup>6</sup> HEK-293 cells were stained with 0.2ug LENG8 antibody (13336-1-AP, red) and control antibody (blue). Cells were fixed with 4% PFA and permeabilized with 0.1% Triton X-100. Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.



Jurkat cells were subjected to SDS PAGE followed by western blot with 13336-1-AP (LENG8 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.