

Nur für Forschungszwecke

# FHIT Polyklonaler Antikörper

Katalog-Nr.:13362-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
13362-1-AP

Größe:  
150ul, Konzentration: 393 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG4181

GenBank-Zugangsnummer:  
BC032336

GeneID (NCBI):  
2272

Vollständiger Name:  
fragile histidine triad gene

Berechnete Masse:  
147 aa, 17 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:  
ELISA

Getestete Reaktivität:  
Human, Maus, Ratte

## Hintergrundinformationen

The candidate tumor suppressor gene, FHIT(Fragile histidine triad protein), encompasses the common human chromosomal fragile site at 3p14.2, the hereditary renal cancer translocation breakpoint, and cancer cell homozygous deletions. Fhit hydrolyzes dinucleotide 5',5'''-P<sub>1</sub>,P<sub>3</sub>-triphosphate in vitro and mutation of a central histidine abolishes hydrolase activity(PMID:9391102). FHIT has a major role in regulating beta-catenin-mediated gene transcription.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

