

Nur für Forschungszwecke

MEI1 Polyklonaler Antikörper

Katalog-Nr.:13456-1-AP



Allgemeine Informationen

Katalog-Nr.:
13456-1-AP

Größe:
150ul, Konzentration: 750 µg/ml von
Nanodrop und 267 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG4253

GenBank-Zugangsnummer:
BC032248

GeneID (NCBI):
150365

Vollständiger Name:
meiosis inhibitor 1

Berechnete Masse:
894 aa, 99 kDa

Beobachtete Masse:
63 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : Maushirngewebe,

Hintergrundinformationen

MEI1 may play a role in meiosis during spermatogenesis. It is required for normal meiotic chromosome synapsis. MEI1 may be involved in the formation of meiotic double-strand breaks (DSBs) in spermatocytes. This antibody recognizes isoform 1, 2, 3, 5, and 7 of MEI1 gene with MW 141 kDa, 67-72 kDa,49-54 kDa and 26 kDa .

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse brain tissue were subjected to SDS PAGE followed by western blot with 13456-1-AP (MEI1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.