

Nur für Forschungszwecke

REG1B Polyklonaler Antikörper

Katalog-Nr.: 13519-1-AP



Allgemeine Informationen

Katalog-Nr.:
13519-1-AP

Größe:
150ul, Konzentration: 227 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG4117

GenBank-Zugangsnummer:
BC027895

GeneID (NCBI):
5968

Vollständiger Name:
regenerating islet-derived 1 beta

Berechnete Masse:
166 aa, 19 kDa

Beobachtete Masse:
19 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



PC-3 cells were subjected to SDS PAGE followed by western blot with 13519-1-AP (REG1B antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.