

Nur für Forschungszwecke

ZNF443 Polyklonaler Antikörper

Katalog-Nr.: 13667-1-AP



Allgemeine Informationen

Katalog-Nr.:

13667-1-AP

Größe:

150ul, Konzentration: 200 µg/ml von Nanodrop und 133 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG4585

GenBank-Zugangsnummer:

BC032753

GeneID (NCBI):

10224

Vollständiger Name:

zinc finger protein 443

Berechnete Masse:

671 aa, 78 kDa

Beobachtete Masse:

72-78 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: Maus-Eierstockgewebe, Mauslungengewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse ovary tissue were subjected to SDS PAGE followed by western blot with 13667-1-AP (ZNF443 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.