

Nur für Forschungszwecke

CA14 Polyklonaler Antikörper

Katalog-Nr.:13736-1-AP



Allgemeine Informationen

Katalog-Nr.:
13736-1-AP

Größe:
150ul, Konzentration: 300 µg/ml von
Nanodrop und 227 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG4673

GenBank-Zugangsnummer:
BC034412

GeneID (NCBI):
23632

Vollständiger Name:
carbonic anhydrase XIV

Berechnete Masse:
337 aa, 38 kDa

Beobachtete Masse:
45-50 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2400

Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : Maus-Augengewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse eye tissue were subjected to SDS PAGE followed by western blot with 13736-1-AP (CA14 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.