

Nur für Forschungszwecke

# JIP2 Polyklonaler Antikörper

Katalog-Nr.:13912-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
13912-1-AP

Größe:  
150ul , Konzentration: 193 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG4902

GenBank-Zugangsnummer:  
BC047527

GeneID (NCBI):  
23542

Vollständiger Name:  
mitogen-activated protein kinase 8  
interacting protein 2

Berechnete Masse:  
824 aa, 88 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:  
IHC,ELISA

Getestete Reaktivität:  
Human

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil  
Lagerungspuffer:  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human testis using 13912-1-AP (JIP2 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis using 13912-1-AP (JIP2 antibody) at dilution of 1:50 (under 40x lens).