

Nur für Forschungszwecke

TSPYL1 Polyklonaler Antikörper

Katalog-Nr.:13932-1-AP



Allgemeine Informationen

Katalog-Nr.:
13932-1-AP

Größe:
150ul, Konzentration: 500 µg/ml von
Nanodrop und 280 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG4961

GenBank-Zugangsnummer:
BC048969

GeneID (NCBI):
7259

Vollständiger Name:
TSPY-like 1

Berechnete Masse:
49 kDa

Beobachtete Masse:
49 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung
Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human, Maus

Positivkontrollen:

WB : Mauslungengewebe, Mausherzgewebe,
Mausnierengewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



mouse lung tissue were subjected to SDS PAGE followed by western blot with 13932-1-AP (TSPYL1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.