

Allgemeine Informationen

Katalog-Nr.:

14019-1-AP

Größe:

150ul, Konzentration: 353 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG5133

GenBank-Zugangsnummer:

BC057837

GeneID (NCBI):

126382

Vollständiger Name:

nuclear receptor 2C2-associated
protein

Berechnete Masse:

16 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

NR2C2AP, also named as NR2C2, TRA16, PGR, PSBPC1, NR3C1, AR and OXT, belongs to the NR2C2AP family. It may act as a repressor of NR2C2-mediated transactivation by suppressing the binding between NR2C2/TR4 and the TR4-response element in target genes. NR2C2AP is mainly localized in the nucleus of cells as a repressor to suppress TR4-mediated transactivation. It can also block the interaction between TR4 and TR4 ligand-binding domain through interacting with TR4-DNA-binding and ligand-binding domains. It may function as a novel repressor to selectively suppress the TR4-mediated transactivation.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

