

Nur für Forschungszwecke

PRR13 Polyklonaler Antikörper

Katalog-Nr.: 14101-1-AP

1 Publikationen



Allgemeine Informationen

Katalog-Nr.:
14101-1-AP

Größe:

150ul, Konzentration: 400 µg/ml von
Nanodrop und 260 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG5280

GenBank-Zugangsnummer:

BC066943

GeneID (NCBI):

54458

Vollständiger Name:

proline rich 13

Berechnete Masse:

15.2 kDa

Beobachtete Masse:

~20 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human

Zitierte Arten:

Human

Positivkontrollen:

WB: K-562-Zellen, MCF-7-Zellen

Hintergrundinformationen

PRR13, also named as TXR1, is reported to be a key regulator of the resistance to cytostatica by decreasing the copy number of the proapoptotic gene thrombospondin-1. PRR13 has some isoforms with MW 10 kDa, 15-17 kDa and 19-23 kDa.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Yun Fu	32606750	Onco Targets Ther	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 14101-1-AP (PRR13 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.