

Nur für Forschungszwecke

IL-18BP Polyklonaler Antikörper

Katalog-Nr.: 14153-1-AP



Allgemeine Informationen

Katalog-Nr.: 14153-1-AP	GenBank-Zugangsnummer: BC044215	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 300 µg/ml von Nanodrop;	GeneID (NCBI): 10068	Empfohlene Verdünnungen: WB 1:500-1:1000
Wirt: Kaninchen	Vollständiger Name: interleukin 18 binding protein	
Isotyp: IgG	Berechnete Masse: 194 aa, 21 kDa	
Immunogen Katalognummer: AG5335	Beobachtete Masse: 21 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : A549-Zellen, HeLa-Zellen
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen

IL18BP, also named as adekinig-alfa, has 3 isoforms (Isoform A 21kd; B 12kd and D 17kd). Isoform A binds to IL-18 and inhibits its activity. IL18BP functions as an inhibitor of the early TH1 cytokine response. Highly glycosylated of IL18BP will migrate to be 40kd. A 58kd protein is IL18BP crosslinked IL18.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 14153-1-AP (IL-18BP antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.