

Nur für Forschungszwecke

# Kv1.2 Polyklonaler Antikörper

Katalog-Nr.: 14235-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
14235-1-AP

Größe:  
150ul, Konzentration: 600 µg/ml von  
Nanodrop und 300 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG5477

GenBank-Zugangsnummer:  
BC043564

GeneID (NCBI):  
3737

Vollständiger Name:  
potassium voltage-gated channel,  
shaker-related subfamily, member 2

Berechnete Masse:  
499 aa, 57 kDa

Beobachtete Masse:  
66-75 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:3000  
IP 0.5-4.0 µg für IP und 1:500-1:1000  
für WB  
IHC 1:50-1:500

## Anwendungen

Geprüfte Anwendungen:

IHC, IP, WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

**Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (\*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.**

Positivkontrollen:

WB: A549-Zellen, HeLa-Zellen, HepG2-Zellen, humanes Herzgewebe, MCF-7-Zellen, Multizellen/-gewebe

IP: HeLa-Zellen,

IHC: Maushirngewebe,

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Ausgewählte Validierungsdaten




Various lysates were subjected to SDS PAGE followed by western blot with 14235-1-AP (Kv1.2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 14235-1-AP (Kv1.2 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 14235-1-AP (Kv1.2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



IP Result of anti-Kv1.2 (IP:14235-1-AP, 4ug; Detection:14235-1-AP 1:500) with HeLa cells lysate 5200ug.