

Nur für Forschungszwecke

FAM175A Polyklonaler Antikörper

Katalog-Nr.:14366-1-AP



Allgemeine Informationen

Katalog-Nr.:
14366-1-AP

Größe:
150ul, Konzentration: 600 µg/ml von
Nanodrop und 333 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG5688

GenBank-Zugangsnummer:
BC039573

GeneID (NCBI):
84142

Vollständiger Name:
family with sequence similarity 175,
member A

Berechnete Masse:
47 kDa

Beobachtete Masse:
46 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:3000

IF 1:20-1:200

Anwendungen

Geprüfte Anwendungen:

IF, WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB : HEK-293T-Zellen, HeLa-Zellen, humanes
Hodengewebe, Jurkat-Zellen, MCF-7-Zellen, Raji-
Zellen

IF : MCF-7-Zellen,

Hintergrundinformationen

The family with sequence similarity 175 member A gene (FAM175A; also known as ABXAS1, CCDC98 and ABRA1), a member of the DNA repair family, contributes to the BRCA1 (BRCA1 DNA repair associated)-dependent DNA damage response and is associated with age at natural menopause.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HEK-293T cells were subjected to SDS PAGE followed by western blot with 14366-1-AP (FAM175A antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed MCF-7 cells using 14366-1-AP (FAM175A antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).