

Nur für Forschungszwecke

TAX1BP1 Polyklonaler Antikörper

Katalog-Nr.:14424-1-AP

Vorgestelltes Produkt

11 Publikationen



Allgemeine Informationen

| | | |
|--|--|--|
| Katalog-Nr.: 14424-1-AP | GenBank-Zugangsnummer: BC050358 | Reinigungsmethode: Antigen-Affinitätsreinigung |
| Größe: 150ul, Konzentration: 1000 µg/ml von8887 Nanodrop; | GeneID (NCBI): Vollständiger Name: Tax1 (human T-cell leukemia virus type I) binding protein 1 | Empfohlene Verdünnungen: WB 1:1000-1:5000 IF 1:10-1:100 |
| Wirt: Kaninchen | Berechnete Masse: 91 kDa | |
| Isotyp: IgG | Beobachtete Masse: 68 kDa, 85 kDa | |
| Immunogen Katalognummer: AG5789 | | |

Anwendungen

| | |
|---|--|
| Geprüfte Anwendungen: IF, WB, ELISA | Positivkontrollen: WB : A549-Zellen, HeLa-Zellen, HepG2-Zellen IF : HepG2-Zellen, |
| In Publikationen genannte Anwendungen: IF, WB | |
| Getestete Reaktivität: Human, Maus | |
| Zitierte Arten: Human, Maus | |

Hintergrundinformationen

TAX1BP1 (also named as T6BP and TXBP151) inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. It is degraded by caspase-3-like family proteins upon TNF-induced apoptosis. TAX1BP1 may also play a role in the pro-inflammatory cytokine IL-1 signaling cascade. TAX1BP1 binds TRAF6 but not TRAFs 1-5. TAX1BP1 has 4 isoforms with MW 65-68 kDa and 86-91 kDa. Western blot analysis showed that TAX1BP1 is expressed as an 86 kDa or a 68 kDa protein.

Bemerkenswerte Veröffentlichungen

| Verfasser | Pubmed ID | Journal | Anwendung |
|--------------|-----------|------------------------|-----------|
| Jing-Ying Li | 36274483 | Int Immunopharmacol | WB |
| Jiajia Pan | 35618171 | Fish Shellfish Immunol | WB |
| Wenxian Wu | 33779513 | Autophagy | WB |

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 14424-1-AP (TAX1BP1 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of HepG2 cells using 14424-1-AP (TAX1BP1 antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).