

Nur für Forschungszwecke

BAK Polyklonaler Antikörper

Katalog-Nr.: 14673-1-AP

34 Publikationen



Allgemeine Informationen

Katalog-Nr.:
14673-1-AP

Größe:
150ul, Konzentration: 300 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG6354

GenBank-Zugangsnummer:
BC004431

GeneID (NCBI):
578

Vollständiger Name:
BCL2-antagonist/killer 1

Berechnete Masse:
23 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

In Publikationen genannte Anwendungen:
IF, IHC, WB

Getestete Reaktivität:
Human, Maus, Ratte

Zitierte Arten:
Human, Maus, Ratte

Hintergrundinformationen

BAK1, also named as BAK, BCL2L7 and CDN1, belongs to the Bcl-2 family. In the presence of an appropriate stimulus, BAK1 accelerates programmed cell death by binding to, and antagonizing the anti-apoptotic action of BCL2 or its adenovirus homolog E1B 19k protein. Low micromolar levels of zinc ions inhibit the promotion of apoptosis. (PMID:17157251) BAK1 gene product primarily enhances apoptotic cell death following an appropriate stimulus. It can inhibit cell death in an Epstein-Barr virus-transformed cell line.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Xufeng Tao	25083618	Transplantation	WB
Yanan Li	27663303	BMC Microbiol	WB
Zuowei Pei	34520392	Aging (Albany NY)	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

