

Nur für Forschungszwecke

C14orf2 Polyklonaler Antikörper

Katalog-Nr.: 14704-1-AP



Allgemeine Informationen

Katalog-Nr.:
14704-1-AP

Größe:

150ul, Konzentration: 400 µg/ml von
Nanodrop und 133 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG6386

GenBank-Zugangsnummer:

BC001944

GeneID (NCBI):

9556

Vollständiger Name:

chromosome 14 open reading frame 2

Berechnete Masse:

7 kDa

Beobachtete Masse:

7 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

IP 0.5-4.0 µg für IP und 1:200-1:1000

für WB

IHC 1:20-1:200

Anwendungen

Geprüfte Anwendungen:

IHC, IP, WB, ELISA

Getestete Reaktivität:

Human

**Hinweis-IHC: Antigendemaskierung mit TE-
Puffer pH 9,0 empfohlen. (*) Wahlweise
kann die Antigendemaskierung auch mit
Citratpuffer pH 6,0 erfolgen.**

Positivkontrollen:

WB: Jurkat-Zellen, humanes Hirngewebe

IP: Jurkat-Zellen,

IHC: humanes Herzgewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.


Ausgewählte Validierungsdaten




Jurkat cells were subjected to SDS PAGE followed by western blot with 14704-1-AP (C14orf2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



IP Result of anti-C14orf2 (IP:14704-1-AP, 4 μ g; Detection:14704-1-AP 1:300) with Jurkat cells lysate 3360 μ g.



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 14704-1-AP (C14orf2 Antibody) at dilution of 1:50 (under 10 \times lens).



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 14704-1-AP (C14orf2 Antibody) at dilution of 1:50 (under 40 \times lens).