

Nur für Forschungszwecke

ANGPTL1 Polyklonaler Antikörper

Katalog-Nr.: 14709-1-AP

1 Publikationen



Allgemeine Informationen

Katalog-Nr.: 14709-1-AP	GenBank-Zugangsnummer: BC050640	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 133 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 9068	
Wirt: Kaninchen	Vollständiger Name: angiopoietin-like 1	
Isotyp: IgG	Berechnete Masse: 55 kDa	
Immunogen Katalognummer: AG6403	Beobachtete Masse: 55 kDa	

Anwendungen

Geprüfte Anwendungen:
ELISA

In Publikationen genannte Anwendungen:
WB

Getestete Reaktivität:
Human

Zitierte Arten:
Hamster

Hintergrundinformationen

Angiopoietin-like protein (Angptl) 1, also known as angiopoietin-related protein 1 (ARP1) and angioarrestin, belongs to the angiopoietin-like family of proteins (Angptls). It is a secretory glycoprotein that shares two basic structural motifs with the angiopoietin family: an N-terminal coiled-coil domain and a C-terminal fibrinogen-like domain. Angptl1 is widely expressed in adult tissues with mRNA levels highest in highly vascularized tissues. It is related to development of the connective tissue and cartilage. This protein was found to be a secretory protein that does not act as an endothelial cell mitogen in vitro.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Akhter Shahina S	23483483	Protein Cell	WB

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

