

WDR4 Polyklonaler Antikörper

Katalog-Nr.: 14766-1-AP

Allgemeine Informationen

Katalog-Nr.:

14766-1-AP

Größe:

150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG6531

GenBank-Zugangsnummer:

BC001074

GeneID (NCBI):

10785

Vollständiger Name:

WD repeat domain 4

Berechnete Masse:

30 kDa

Beobachtete Masse:

30 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

A group of functionally diverse regulatory proteins in eukaryotes, defined by a conserved core of 4 or more modular repeat units, is known as the WD repeat family. Proteins belonging to the WD repeat family are involved in a number of cellular processes, including cell cycle progression, apoptosis, signal transduction, and gene regulation. WD repeat-containing protein 4 (WDR4) belongs to the WD repeat family, and is required for the formation of N(7)-methylguanine at position 46 (m7G46) in tRNA. It acts to stabilize and induce conformation change of the catalytic subunit in the complex. There exists three isoforms, while isoform 1, 2 are 45kDa band, and isoform 3 is 30 kDa band.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

