

Nur für Forschungszwecke

IL-17D Polyklonaler Antikörper

Katalog-Nr.: 14934-1-AP



Allgemeine Informationen

Katalog-Nr.:
14934-1-AP

Größe:
150ul, Konzentration: 253 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG6765

GenBank-Zugangsnummer:
BC036243

GeneID (NCBI):
53342

Vollständiger Name:
interleukin 17D

Berechnete Masse:
22 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

IL17D is a secreted cytokine with homology to the IL-17 family of proteins. IL17D is preferentially expressed in skeletal muscle, brain, adipose tissue, heart, lung, and pancreas. The treatment of endothelial cells with IL17D has been shown to stimulate the production of other cytokines including IL6, IL8 and CSF2/ GM-CSF. The increased expression of IL8 induced by IL17D was found to be NF-kappa B-dependent.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

