

SNX11 Polyklonaler Antikörper

Katalog-Nr.: 14953-1-AP

Allgemeine Informationen

Katalog-Nr.:

14953-1-AP

Größe:

150ul, Konzentration: 600 µg/ml von Nanodrop und 287 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG6791

GenBank-Zugangsnummer:

BC000768

GeneID (NCBI):

29916

Vollständiger Name:

sorting nexin 11

Berechnete Masse:

30 kDa

Beobachtete Masse:

35 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB: HeLa-Zellen,

Hintergrundinformationen

SNX11 belongs to the sorting nexin family of proteins which contain a phox (PX) domain. PX domain is a phosphoinositide binding domain. SNX family members are involved in intracellular trafficking. SNX11 protein lacks C-terminal coiled-coil motifs found in several other SNX family members.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HeLa cells were subjected to SDS PAGE followed by western blot with 14953-1-AP (SNX11 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.