

Nur für Forschungszwecke

# RIMS3 Polyklonaler Antikörper

Katalog-Nr.:15088-1-AP



## Allgemeine Informationen

**Katalog-Nr.:**

15088-1-AP

**Größe:**

150ul, Konzentration: 200 µg/ml von Nanodrop und 173 µg/ml durch die Bradford-Methode mit BSA als Standard;

**Wirt:**

Kaninchen

**Isotyp:**

IgG

**Immunogen Katalognummer:**

AG7171

**GenBank-Zugangsnummer:**

BC003103

**GeneID (NCBI):**

9783

**Vollständiger Name:**

regulating synaptic membrane exocytosis 3

**Berechnete Masse:**

33 kDa

**Beobachtete Masse:**

40 kDa

**Reinigungsmethode:**

Antigen-Affinitätsreinigung

**Empfohlene Verdünnungen:**

WB 1:500-1:1000

## Anwendungen

**Geprüfte Anwendungen:**

WB, ELISA

**Getestete Reaktivität:**

Human, Maus, Ratte

**Positivkontrollen:**

WB: HL-60-Zellen, Maushirngewebe

## Hintergrundinformationen

### Lagerung

**Lagerungsbedingungen:**

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

**Lagerungspuffer:**

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



HL-60 cells were subjected to SDS PAGE followed by western blot with 15088-1-AP (RIMS3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.