

Nur für Forschungszwecke

DLX2 Polyklonaler Antikörper

Katalog-Nr.:15174-1-AP



Allgemeine Informationen

Katalog-Nr.:

15174-1-AP

Größe:

150ul , Konzentration: 800 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG4900

GenBank-Zugangsnummer:

BC032558

GeneID (NCBI):

1746

Vollständiger Name:

distal-less homeobox 2

Berechnete Masse:

34 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

DLX2 belongs to the DLX family of homeobox transcription factors (Dlx1, Dlx2, Dlx5 and Dlx6), which is expressed in overlapping domains at different stages of cell differentiation in the subpallium and controls differentiation of GABAergic neurons. DLX2 may play a role regulation in the development of the ventral forebrain. Also it may involve in craniofacial patterning and morphogenesis. DLX2 can form homodimer with DLX5(-34kd)

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

